

Paretor – Predictable Delivery Protocols & Formulations

Mineral Nanoparticle Colloid | Cell transfection | Cargo Carrier

For research use only. Not for use in Diagnostic nor Therapeutic Procedures.

Up to 65%
transfection
efficiency in
available cells
in 2 hours*

APPLICATIONS:

Cytosolic delivery of cargo including siRNA, microRNA, ASO, Antagomirs up to 30nt, small molecules, proteins up to 15KDa

STABILITY:

One month at 4 °C, 21 °C.
Nine months at -20 °C
Freeze-thaw stable, thaw on ice

Size:	40-45nm
PDI:	<0.3
Control Charge:	-27mV
Carrier Charge:	-15mV
Polymer:	Dextran
Cargo Conjugation:	Thiol-reactive
Surface Logic:	Lipid Head
Surface Passivation:	Zwitterionic
Magnetic Saturation:	< 1emu/g

Particles Received:

PA-001 base construct (10-001, 10-002, 10-003, 10-004, 10-005, and 10-005)

Derived from PA-001

Concentration: **1.0 × 10¹¹ particles/mL**

Total volume: **1.0 mL** (4 × 250 µL aliquots)

Storage Buffer System:

0.5 mM Sodium Gluconate

2 mM Sodium Citrate

5% Trehalose

pH neutral, isotonic, RNase-free

Surface Charge:

Control (10-001): -27 mV

Carrier (10-002): -15 mV

* Some cell types exhibit **charge sensitivity**; expect differences in uptake depending on metabolic state and confluency.

Storage & Stability:

Freeze-thaw stable

Stability at -20 °C: validated to nine months

Stability at 4 °C: Validated to 4 weeks

Store dark + dry regardless of temperature

Thaw ON TOP of ice (not submerged) to prevent uncontrolled warming

Inspect for particulates (should be none). If present, discard vial and notify.

Functional Chemistry:

Carrier particles include maleimide group adducts for thiol-reactive surface conjugation

Zwitterionic surface logic; ligand shell is stable and non-aggregating

Avoid moisture and light to preserve functional groups

* Demonstrated in BT20 cells. Your results may vary. Not approved for clinical use. User responsible for safety, compliance.